



# mosquito® med Head Lice Shampoo10

## Instructions for use

Before using mosquito® med Head Lice Shampoo10, please carefully read the following instructions for use. If you have any questions, please consult your doctor or chemist.

## Ingredients

Mineral oil, foaming agent, perfume

## Areas of application

mosquito® med Head Lice Shampoo10 is used for the physical treatment of head hair infested with head lice. If it is used correctly, it eliminates head lice and nits.


## Physical active principles

The combination of active ingredients in mosquito® med Head Lice Shampoo10 blocks the airways (trachea) of head lice and nits so that they suffocate. It also dissolves the protective wax layer on the louse's carapace, causing it to dry out.

## Counter-indicators

- Do not use mosquito® med Head Lice Shampoo10 if you are aware of hypersensitivity or allergies to mineral oils. If you encounter skin eruptions or redness, stop treatment immediately. Do not use the product if the scalp is irritated or damaged.
- mosquito® med Head Lice Shampoo10 is not suitable for children aged less than 12 months. Due to lack of empirical evidence, we advise against using the product when pregnant or while breastfeeding.

## ⚠ Warnings

- Keep the product out of reach of children.
- mosquito® med Head Lice Shampoo10 is intended for external use only.
- Avoid contact with the eyes and mucous membranes. In rare cases, mosquito® med Head Lice Shampoo10 can cause serious irritation to the eyes.
- If shampoo should get into the eyes, thoroughly rinse them **immediately** with water.
- In the event of accidental ingestion, immediately seek medical attention. 
- Because some of the ingredients are flammable, persons whose hair has been treated should not stand too close to naked flames. Do not smoke during use.
- It is impossible to exclude the possibility that mosquito® med Head Lice Shampoo10 may react with various hair care products and hair colorants.
- Do not use mosquito® med Head Lice Shampoo10 if you are aware of hypersensitivity or allergies to mineral oils.
- If you accidentally spill mosquito® med Head Lice Shampoo10, it can cause extremely slippery conditions. It is best to wash the affected area directly.

**mosquito® med Head Lice Shampoo10 is applied to dry hair**

## Storage guidelines

Do not store at over 25 degrees. Do not use mosquito® med Head Lice Shampoo10 after the expiration date.



Head Lice  
Shampoo10  
100 ml  
PZN 10415469

Head Lice  
Shampoo10  
200 ml  
PZN 10415475

## Treat the hair in three steps:

### Step 1: Apply the shampoo to dry hair

The amount to be applied largely depends on the length and thickness of the hair. Spread the shampoo in the dry hair and make sure that all the hair along its complete length is completely saturated with the shampoo. Likewise, all areas of the scalp must be treated. Take particular care to ensure that the areas behind the ears, at the nape of the neck and the temples are treated with sufficient amounts of shampoo.

**mosquito® med Head Lice Shampoo10 can only work in areas where it is applied.** Make sure that no shampoo gets in the eyes, for example by covering them with a wash cloth. Leave the shampoo in the hair for at least 10 minutes without covering it (with a towel, hood, etc.). This is very important for a successful treatment.



### Step 2: Wash out the shampoo

Apply a small amount of water to the hair and massage the shampoo into a lather. Then rise the foam out using luke warm water. Make sure not to get the foam in the eyes. Rinse the hair thoroughly until no shampoo remains. There is no need to treat the hair with a normal commercial shampoo afterwards.

**Tip:** Rinsing increases the combability of the hair. (See step 3)

### Step 3: Comb the hair using a lice comb

A head lice comb is a special type of comb with tightly packed teeth. Start by applying the comb directly at the scalp and comb out the hair strand-by-strand. Every time you comb the hair, clean head lice and nits from the comb with an absorbent cotton cloth or kitchen towel. Make sure that the hair is damp while you comb it out.

**Repeat the application (steps 1 , 2 and 3) after 7–10 days.**

### What further steps should you take?

- Notify the school, crèche or contact person that there has been a head lice outbreak.
- As a precautionary measure, treat all family members with mosquito® med Head Lice Shampoo10.
- Wash all clothing, bed linen and towels used at a minimum of 60 degrees or use mosquito® Detergent, which can be used at 30 degrees or more. Non-washable objects can be treated with the mosquito® Fabric Spray or can be stored in a tightly sealed plastic bag for at least 4 days.

### Disposal instructions

Send empty containers for recycling.

Information date: July 2015

Distributed by:



WEPA Apothekenbedarf GmbH & Co KG  
56204 Hillscheid  
www.wepa-dieapothekenmarke.de



Head Lice Shampoo10 manufacturer:  
Oystershell NV  
Booiebos 24  
B-9031 Drongen (Belgium)



07020363  
1086/020097/15D

[www.mosquito-laese.de](http://www.mosquito-laese.de)