




How to apply mosquito® Head Lice 2in1 Fluid

Activating the bottle



1. Hold the bottle firmly in your hand and place your thumb on point 1.
2. Push the cap forward. The cap is now positioned at an angle. The active ingredient liquid is released.
3. Shake the bottle thoroughly for 10 seconds.
4. Break off the tip of the cap. The product is now ready for use. Don't throw the tip away. You must close the bottle again afterwards.

Apply to dry hair!

Step 1: Spread the fluid on the dry hair

Put a towel around the person's neck and shoulders during treatment. Apply a generous amount of the fluid to the dry hair. The hair and scalp must be completely saturated with the fluid. Massage the product into hair and scalp thoroughly, paying particular attention to the areas behind the ears and at the nape of the neck. The fluid must now act for 15 minutes (no longer).



Step 2: Comb the hair using a lice comb

During this 15-minute application time the hair can be easily combed. We recommend combing out the hair with a lice comb during the 15-minute application time. If necessary, comb the hair with a normal brush to untangle it. Then proceed as follows: Separate the hair into four sections. Comb each section with the lice comb strand by strand. Start as close as possible to the hairline and comb it all the way to the hair-ends. In this way you remove the lice/nits. Any strand that is free from lice/nits can be put aside. Remove all lice/nits from the comb with a cotton pad, handkerchief or similar. Alternatively, you can also comb out the hair after it has been rinsed.




Step 3: Rinsing the hair

After the 15-minute application time, rinse out the fluid. To do this, wash the hair thoroughly with a standard shampoo and rinse well.



Step 4: Follow-up check

For every head lice agent the Robert Koch Institute recommends a repeat treatment after 9 days (+/- 24 hours). When used as intended, mosquito® Head Lice Fluid 2in1 is already effective after a single treatment. However, a follow-up check should be carried out 9 days (+/- 24 hours) after the treatment. If you discover live head lice, you should repeat the treatment as described above.



Shake the bottle before each use.



Head Lice 2in1 Fluid
100 ml
PZN 16362998
Also available in 200 ml
family size
PZN 16363006



WEPA Apothekenbedarf GmbH & Co KG
D-56204 Hillscheid
www.wepa-dieapothekenmarke.de

You will also find a product application video at www.mosquito-laese.de under „Service“ or here:

